



GOLDROACH

PROTECTS FROM THE INSIDE OUT



REUSABLE AND WASHABLE ELEMENT HIDER

- Fabric 5168-neo
- BLOOD BARRIER treatment
- C6 TecHnology
- APPLIFRESH BC with static permanent microbial agent.
- AITEX APPROVED FABRIC WITH OEKO-TEX 100 CERTIFICATE

SIZES (CM)

SMALL	50-54
MEDIUM	55-59
LARGE	60-64

From tip of nose to tip of nose passing through the nape

REUSABLE AND WASHABLE ELEMENT HIDER

5168-NEO fabric formed by two layers assembled in its bond, with a weight of 280gr / m² and a composition of 92% PES, 8% EA.

BLOOD BARRIER formulations:

- **APPLIGUARD DPS:** Composite Fluorocarbon based on C6 technology for fluid repellency.
- **APPLIFRESH BC:** Static permanent microbial agent (does not migrate).
- **RETICULANTES BI-25:** Helps the fluorocarbon stay longer on the fiber.
- **ANTIESTATICO AF:** Eliminates static electricity and thus prevents the adhesion of particles.

THE FORMULATIONS FOR OBTAINING FLUID-REPELLENT AND ANTIBACTERIAL FINISHES (“BLOOD BARRIER”) BEGIN TO DETERIORATE FROM 25 WASHES AT 40OC (HL40OC) OR 5 WASHES AT 60OC (HL60OC)

LEGAL DISCLAIMER

“WARNING: THIS DEVICE IS NOT A MEDICAL DEVICE WITHIN THE MEANING OF DIRECTIVE 93/42 OR REGULATION EU / 2017/745, NOR IS IT A PERSONAL PROTECTIVE EQUIPMENT WITHIN THE MEANING OF REGULATION EU / 2016/425”;